



Decavita

**ALL YOU
NEED FOR
ALIGNERS**

GC

DentaSonic
Swiss Proximal Restorative Clean-up Tools

TASK

oscident[®]

EKULF

NOTES

GC BONDING MATERIALS



GC ALIGNER CONNECT 90-4100-0000

GC Aligner Connect was specifically developed for orthodontic aligner attachments. GC Aligner Connect features enhanced aesthetics that imitate enamel color shades and high durability which limits long term wear or potential detachment during treatment.



GC ORTHO CONNECT FLOW 90-3100-0000

Light cure orthodontic adhesive for lingual retainers, stops and aligner attachments – easy to handle. GC Ortho Connect Flow's primer is integrated into the paste, allowing to save valuable time in the bonding process.













GC ORTHO CONNECT 90-1100-0000




A new orthodontic light-cured adhesive for brackets and bondable buccal tubes. GC Ortho Connect's primer is integrated into the paste, allowing to save valuable time in the bonding process.



BONDABLE LINGUAL ATTACHMENTS

Bondable lingual buttons are used in dental movements such as retractions, rotations, pull and anchoring systems. This retains the elastic, wires, chains and ligatures.

| Bondable Buttons | | P/N | |
|---|---|----------------------------------|--------------|
|  |  | Flat base, pouch of 10 | 42-0020-0000 |
| |  | Curved base, pouch of 10 | 42-0010-0000 |
|  |  | Curved base Rhodium, pouch of 10 | 42-0010-RH00 |

| Bondable Cleats | | P/N | |
|---|---|--------------------------|--------------|
|  |  | Flat base, pouch of 10 | 42-0060-0000 |
| |  | Curved base, pouch of 10 | 42-0050-0000 |
|  |  | Eyelet, pouch of 10 | 42-0070-0000 |

| Surgical buttons | | P/N | |
|---|---|--------------------------|--------------|
|  |  | Curved base, pouch of 10 | 42-0030-0000 |
|  |  | Flat base, pouch of 10 | 42-0040-0000 |

| Gold Chain by 5 | | P/N | |
|---|------------------------------|--------------|--|
|  | Round base, pouch of 5 | 42-3150-0000 | |
|  | Rectangular base, pouch of 5 | 42-3160-0000 | |



METAL CHAIN WITH
ROUND BASE (2 psc/pack)
ILB023-11



CERAMIC BUTTON FREDDY (10 psc/pack)
Curved Base
44-1010-0000

The Freddy ceramic button is the perfect clear utility button for all aesthetic appliances.

- Works lingually as well as labially
- Mesh structure of base makes a strong bonding
- Bond well with most common adhesives
- Anatomy and size base fits all teeth
- Enough space around tie wing for elastics
- Maintains the transparent look with aligners as anchorage for elastics



CERAMIC HOOKS (10 psc/pack)
Curved Base
FCB-00-01HK



PLASTIC BUTTONS (10 psc/pack)
Curved Base
1545



ALIGNER CHEWIES BITE ROLLS


1560 – Tasteless

Chewies help the patient seat the aligner into place securely and comfortably. They also stimulate the periodontal ligaments to aid effective tooth movement. Each package contains 10 bags of 2 chewies.



ORCHESTRA INTRA AND EXTRAORAL ELASTICS

Our range of intra and extraoral elastics are made from the highest quality medical grade latex. The stated force corresponds to the force three times its original listed diameter. Intraoral elastics are packaged by 5000 pcs. Extraoral elastics are packaged by 2500 pcs. Each bag contain 1 insertion plastic instrument. Latex free elastics are listed on the chart below in green.



JAPANESE ELASTOMERIC CHAIN – GC ORTHO CHAIN

Made from high-density and high-purity materials, this chain is precision-cut with identical intervals between links. GC Ortho Chain exerts light and continuous force throughout its action, which is essential for optimal tooth movement. Available in 2 m spools in grey or clear.

| Color | Short | Middle | Long |
|-------------|--------------|--------------|--------------|
| Transparent | 86-1010-0000 | 86-1020-0000 | 86-1030-0000 |
| Grey | 86-2010-0000 | 86-2020-0000 | 86-2030-0000 |


| EXTRA-ORAL | | | | |
|------------------------------|----------------------------|---|---|-----------------------------|
| Force Size | 2.5 oz (70.88 g) | 4.5 oz (113.40 g) | 6.5 oz (170.10 g) | 8.5 oz (240.98 g) |
| 3 mm (1/8") | 85-0321-0000 Harp | 85-0331-0000 Xylophone | 85-0341-0000 Tuba | |
| 4 mm (3/16") | 85-0421-0000 Tambourine | 85-0431-0000 Piano | 85-0441-0000 Trumpet | 85-0461-0000 Trombone |
| 4 mm (3/16") (sans latex) | | 85-0430-0000 Jazz (latex free) | 85-0440-0000 Rock-n-Roll (latex free) | |
| 6 mm (1/4") | 85-0621-0000 Oboe | 85-0631-0000 Flute | 85-0641-0000 Baritone | 85-0661-0000 French Horn |
| 6 mm (1/4") (sans latex) | | 85-0630-0000 Hip-Hop (latex free) | 85-0640-0000 Classical (latex free) | |
| 8 mm (5/16") | 85-0821-0000 Cello | 85-0831-0000 Clarinet | 85-0841-0000 Saxophone | 85-0861-0000 Guitar |
| 10 mm (3/8") | 85-1021-0000 Cymbals | 85-1031-0000 Tympani | | 85-1061-0000 Violin |

STRIPPING TOOLS

An innovative Interproximal Reduction system that makes interproximal enamel removal/ stripping safer and more accurate without creating sharp corners or subgingival ledges.

Designed for use with all orthodontic treatments, these strips curve and fly along the natural contours of the teeth to avoid cutting lips and gums, prevent excessive enamel removal and to maximize patient comfort.

OSCIDENT IPR STRIPS



| Article | Description | REF-No. |
|-----------------|--------------------------------|---------|
| ● Strip (15 µm) | Polisher, extra fine | 1015 |
| ● Strip (25 µm) | Strip perforated, fine | 1025 |
| ● Strip (40 µm) | Strip perforated, medium | 1040 |
| ● Strip (60 µm) | Strip perforated, coarse | 1060 |
| ● Strip (90 µm) | Strip perforated, extra coarse | 1090 |
| ● Strip (15 µm) | Opener extra fine | 1115 |
| ● Strip (40 µm) | Opener medium | 1140 |



STARTERSETSTRIPS (6) INCL. MEASURING GAUGE

1006

DENTASONIC IPR STRIPS

Refill

2 sides

| Name: | DiaStrip | DiaStrip | DiaStrip | DiaStrip |
|--------------------|---------------------|---------------------|---------------------|---------------------|
| P/N 2 sides, 3pcs: | 55-FDM0-3215 | 55-FDM0-3225 | 55-FDM0-3240 | 55-FDM0-3250 |
| Grain µm: | 15 µm | 25 µm | 40 µm | 60 µm |
| Thickness 2 sides: | 0.10mm | 0.15mm | 0.20mm | 0.30mm |
| Packaging refills: | 3 pieces | 3 pieces | 3 pieces | 3 pieces |



HANDLE DOUBLE STRAIGHT AND 45° ANGLED

1501



DS MEASURING GAUGE

55-G000-1000

Made from Stainless Steel – 0.1, 0.15, 0.2, 0.3, 0.4, 0.5. Measure the exact amount of interproximal reduction achieved.

INTERPROXIMAL REDUCTION PRODUCTS

DENTASONIC IPR KIT



Silent hand piece,
no vibration, no noise.
Non water cooled

IPR HANDPIECE

| | | |
|--|------|----------------------|
| EVA contra angle hand piece Sirona T1 Line 11L | 3334 | with spray and light |
| EVA contra angle hand piece Sirona T3 Line 11 | 3333 | with spray |
| EVA contra angle hand piece NSK TI-Max XSSL | 1210 | with spray and light |
| EVA contra angle hand piece NSK TI-Max X55 | 1211 | with spray |
| W & H WG-69 LT | 1252 | with spray and light |
| W & H WG-69 A | 1253 | with spray |

IPR BURS

Yellow (0.1-0.25mm)
Blue (0.3-0.5mm)





**OSCIDENT ALIGNER
PLIER EMBOSS**

**Retention
1674**

Creates a circular indentation. 1.4mm in diameter.
For increased retention in the retainer and aligner system.



**TASK ALIGNER PLIER
VERTICAL HOLE**

**Rotation
60-733**

Creates a vertical indentation. 1.5mm x 0.75 mm.
For individual tooth rotation.



**TASK ALIGNER PLIER
HORIZONTAL HOLE**

**Torque
60-732**

Creates a horizontal indentation. 0.75mm x 1.5 mm.
For individual root torque.



**TASK ALIGNER PLIER
TEAR DROP HOLE**

60-734

Creates a hook in the gingival margin for the use of elastomerics.



**TASK ALIGNER PLIER
HOLE PUNCH**

60-731

Creates a circular cutout to allow bondable attachments to be used.

ALIGNER PLIER INSTRUMENT SET

1676

The popular choice for premium instrumentation of orthodontists globally. This kit includes all 5 of our Aligner Pliers.

ALIGNER INSTRUMENTS

**ALIGNER REMOVAL
TOOL (30 pcs/ pack)**



MATERIALS FOR ALIGNERSTS

SMILE RITE THERMOFORMING PLASTIC



All Smile Rite Thermoforming plastic sheets, regardless of thickness, are pre-dried and packaged in aluminium bags of 25 pieces. The water absorption by the material has been eliminated as well as the bubble problem during thermoforming has been resolved.

| Hard | Round 120 mm | Round 125 mm | Square 125 mm x 125 mm |
|-----------------|--------------|--------------|------------------------|
| 1,5 mm (x 50) | — | 92-1015-R125 | 92-1015-S125 |
| 1,0 mm (x 100) | 92-1010-R120 | 92-1010-R125 | 92-1010-S125 |
| 0,80 mm (x 100) | — | 92-1008-R125 | — |
| Soft | Round 120 mm | Round 125 mm | Square 125 mm x 125 mm |
| 1,0 mm (x 100) | 92-1110-R120 | 92-1110-R125 | 92-1110-S125 |

SCISSORS FOR ALIGNER FOIL



GOSPET FOIL PREPARATION SET



GOSPET PLASTIC



One appliance
for 90% cases

MULTI FUNCTIONAL APPLIANCE FOR CL II and CL III CORRECTION

MULTI FA regular, short, long
MULTI FA clear
MULTI FA molars



Improvement of
the biomechanical control



Double hinge for an
anatomical positioning
without compromise

Fix the hinge upon
treatment plan
requirements

RETAINER WIRE (10 pcs/pack)
Pre-bended
only 0.17 mm thickness



RELiance BOND-A-BRAID
LINGUAL RETAINER WIRE
(10 pcs/ pack)



TOOLS FOR PATIENT HYGIENE



EKULF ORTHO
IMPLANT KIT BLISTER



EKULF TRAVEL KIT



EKULF TWICE,
JUNIOR/ ADULT BLISTER



EKULF PH SUPREME MIXED,
7 PCS CYLINDRIC, BLISTER



EKULF TWICE
ADVANCED, BLISTER



EKULF PURE SHINE
TOOTHPASTE 75 ML



EKULF MONOBRUSH
SOFT, BLISTER



UAB Decavita
Raminta Rusakevičiūtė
+37069964256
info@decavita.eu